


# Vitapex®

ROOT CANAL MEDICAMENT (Radiopaque)



|                          |              |   |
|--------------------------|--------------|---|
| <b>Calcium hydroxide</b> | <b>30%</b>   | Stimulates “blast” (osteoblast/cementoblast) cells aiding apexogenesis. High PH neutralizes endotoxins produced by anaerobic bacteria |
| <b>Iodoform</b>          | <b>40.4%</b> | Bacteriostatic  |
| <b>Silicone Oil</b>      | <b>22.4%</b> | Increased radiopacity   |
| <b>Inert</b>             | <b>6.9%</b>  | Lubricant, ensures complete coating of canal walls<br>Never hardens   |

## Packaging Storage

2g of paste per syringe + 20 disposable tips  
Store in a dark environment

**Vitapex** can be used as a temporary or permanent root canal filling material. It is ideal for the treatment of infected root canals and for vital pulpotomies in deciduous teeth.

## References

Thomas A.M., et al.

Cohen S., Burns R.  
Eda S., et al.

Nurko C., Garcia-Godoy F.

Elimination of infection in pulpectomized teeth: a short-term study using iodoform paste. J. Endod. 1994 20(5): 233-5  
Pathways of the Pulp. 7th Ed. 1998, Mosby  
Clinical pathological studies on the healing of periapical tissues in aged patients by root canal filling using pastes of Calcium hydroxide added iodoform. Geriodontics 1985 1: 98-104  
Evaluation of a Calcium hydroxide / iodoform paste (V itapex) in root canal therapy for primary teeth.  
J. Clin. Pediatr. Dent 1999 23(4): 289-294




**NEO DENTAL INTERNATIONAL INC.**  
2505 South 320 Street, Suite 625  
Federal Way, WA 98003  
USA

Tel (253) 946-0814 **Toll Free: 1-800-720-0096**  
Fax (253) 946-0885 E-mail: [ndi@neodental-intl.com](mailto:ndi@neodental-intl.com)

**Web: [www.neodental-intl.com](http://www.neodental-intl.com)**

# Vitapex®

ROOT CANAL MEDICAMENT (Radiopaque)



1 prefilled syringe  
containing 2.0g paste  
20 disposable tips #40

## Premixed Calcium Hydroxide & Iodoform Paste for Root Canal Treatment



Vitapex  
Tips Set  
2g Syringe + 20 tips



Disposable Tips  
40 tips/box



Disposable  
Mini Syringes  
20 Syringes & Plungers



Vitapex  
Refill  
2g Syringe

## References

Thomas A.M., et al.

Cohen S., Burns R.  
Eda S., et al.

Nurko C., Garcia-Godoy F.

Elimination of infection in pulpectomized teeth: a short-term study using iodoform paste. J. Endod. 1994 20(5): 233-5  
Pathways of the Pulp. 7th Ed. 1998, Mosby  
Clinical pathological studies on the healing of periapical tissues in aged patients by root canal filling using pastes of Calcium hydroxide added iodoform. Geriodontics 1985 1: 98-104  
Evaluation of a Calcium hydroxide / iodoform paste (V itapex) in root canal therapy for primary teeth.  
J. Clin. Pediatr. Dent 1999 23(4): 289-294

## Indications

- Intracanal medicament
- Apexification
- Exudation control
- Periapical lesions
- Root resorption
- Temporary root canal filling material
- Perforations
- Under developed pulpless teeth

## Benefits

- Helps to quiet “Hot Abscess”
- Disinfects canal
- Antibacterial
- Bacteriostatic
- Radiopaque
- Easy to use
- Promotes apexification
- Aids in calming infection



**NEO DENTAL INTERNATIONAL INC.**

2505 South 320 Street, Suite 625  
Federal Way, WA 98003 USA

**Web: [www.neodental-intl.com](http://www.neodental-intl.com)**

Tel: (253) 946-0814 **Toll Free: 1-800-720-0096**  
Fax: (253) 946-0885 E-mail: [ndi@neodental-intl.com](mailto:ndi@neodental-intl.com)