


Dentalis®

ROOT CANAL SEALER



Calcium Hydroxide
& Iodoform Sealer

Indications

- Root canal sealer
- Antibacterial dressing



Benefits

- Excellent radiopacity
- Antibacterial
- Ideal flow characteristics
- Adheres to root canal wall and obturating points
- Ensures closing of the radicular pulp cavity
- Sets quickly in the mouth

Technical Data

| | |
|---------------------|-----------------------|
| Radiopacity | 4.4mmAL |
| Flow | 31.0mm |
| Solubility | 1.75% |
| Setting time | 5 - 10 minutes |
| Working time | 4 - 6 hours |

Package contains:
powder

Calcium Hydroxide
Iodoform
Zinc Oxide
Others

Liquid

Eugenol
Others

Spoon and Mixing pad

Exceeds International Standards

- **RADIOPAQUE**
- **MEDICATED**
- **OUTSTANDING FLOW**



NEO DENTAL INTERNATIONAL INC.

2505 South 320 Street, Suite 625
Federal Way, WA 98003 USA

Web: www.neodental-intl.com

Tel (253) 946-0814
Fax(253) 946-0885

Toll Free: 1-800-720-0096
E-mail: ndi@neodental-intl.com